

For Research Use Only

# CMC1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60513-3-PBS



## Basic Information

Catalog Number:

60513-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21231

GenBank Accession Number:

BC052644

GeneID (NCBI):

152100

UNIPROT ID:

Q7Z7K0

Full Name:

COX assembly mitochondrial protein  
homolog (S. cerevisiae)

Calculated MW:

106 aa, 12 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B6E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50717-2, CMC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60513-1-PBS. Detection antibody: 60513-3-PBS. Standard:Ag21231. Range: 0.391-100 ng/mL.