

For Research Use Only

CWC15 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60514-3-PBS



Basic Information

Catalog Number:

60514-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG18323

GenBank Accession Number:

BC040946

GeneID (NCBI):

51503

UNIPROT ID:

Q9P013

Full Name:

CWC15 spliceosome-associated protein homolog (S. cerevisiae)

Calculated MW:

229 aa, 27 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A4E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50720-2, CWC15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60514-3-PBS. Detection antibody: 60514-2-PBS. Standard: Ag18323. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50720-4, CWC15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60514-3-PBS. Detection antibody: 60514-4-PBS. Standard: Ag18323. Range: 3.125-100 ng/mL.