

For Research Use Only

BCL7B Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60517-3-PBS



Basic Information

Catalog Number: 60517-3-PBS	GenBank Accession Number: BC000956	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9275	CloneNo.: 2B7A4
Source: Mouse	UNIPROT ID: Q9BQE9	
Isotype: IgG1	Full Name: B-cell CLL/lymphoma 7B	
Immunogen Catalog Number: AG22892	Calculated MW: 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50729-2, BCL7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60517-3-PBS. Detection antibody: 60517-2-PBS. Standard: Ag22892. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50729-3, BCL7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60517-4-PBS. Detection antibody: 60517-3-PBS. Standard: Ag22892. Range: 0.195-100 ng/mL.