

For Research Use Only

AMFR Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60519-4-PBS



Basic Information

Catalog Number: 60519-4-PBS	GenBank Accession Number: BC069197	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 267	CloneNo.: 3A12C7
Source: Mouse	UNIPROT ID: Q9UKV5	
Isotype: IgG1	Full Name: autocrine motility factor receptor	
Immunogen Catalog Number: AG10308	Calculated MW: 643 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50731-3, AMFR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60519-1-PBS. Detection antibody: 60519-4-PBS. Standard: Ag10308. Range: 0.195-50 ng/mL