

For Research Use Only

NTF3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60520-4-PBS



Basic Information

Catalog Number: 60520-4-PBS	GenBank Accession Number: BC069773	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4908	CloneNo.: 2C1H3
Source: Mouse	UNIPROT ID: P20783	
Isotype: IgG2a	Full Name: neurotrophin 3	
Immunogen Catalog Number: AG35836	Calculated MW: 257 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50760-3, NTF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60520-4-PBS. Detection antibody: 60520-1-PBS. Standard: Ag12790. Range: 1.563-200 ng/mL.