

For Research Use Only

# RPS26 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60521-2-PBS



## Basic Information

Catalog Number: 60521-2-PBS	GenBank Accession Number: BC002604	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6231	CloneNo.: 1C2A3
Source: Mouse	UNIPROT ID: P62854	
Isotype: IgG1	Full Name: ribosomal protein S26	
Immunogen Catalog Number: AG6718	Calculated MW: 13 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP5073-1, RPS26 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60521-1-PBS. Detection antibody: 60521-2-PBS. Standard:Ag6718. Range: 0.781-100 ng/mL.