

For Research Use Only

SGPL1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60523-2-PBS



Basic Information

Catalog Number: 60523-2-PBS	GenBank Accession Number: NM_003901	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8879	CloneNo.: 4F12C12
Source: Mouse	UNIPROT ID: O95470	
Isotype: IgG2a	Full Name: sphingosine-1-phosphate lyase 1	
Immunogen Catalog Number: AG26959	Calculated MW: 64 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50739-1, SGLT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60523-1-PBS. Detection antibody: 60523-2-PBS. Standard: Ag26959. Range: 1.563-200 ng/mL.