

For Research Use Only

# CD141/Thrombomodulin Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60524-1-PBS



## Basic Information

Catalog Number:

60524-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_000361.3

GeneID (NCBI):

7056

UNIPROT ID:

P07204

Full Name:

thrombomodulin

Calculated MW:

60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G11C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50741-1, CD141/Thrombomodulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60524-1-PBS. Detection antibody: 60524-2-PBS. Standard:Eg0944. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP50741-2, CD141/Thrombomodulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60524-1-PBS. Detection antibody: 60524-3-PBS. Standard:Eg0944. Range: 0.098-50 ng/mL