

For Research Use Only

CD141/Thrombomodulin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60524-3-PBS



Basic Information

Catalog Number:

60524-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

NM_000361.3

GeneID (NCBI):

7056

UNIPROT ID:

P07204

Full Name:

thrombomodulin

Calculated MW:

60 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3C5D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50741-2, CD141/Thrombomodulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60524-1-PBS. Detection antibody: 60524-3-PBS. Standard:Eg0944. Range: 0.098-50 ng/mL