

For Research Use Only

OS9 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60527-1-PBS



Basic Information

Catalog Number: 60527-1-PBS	GenBank Accession Number: BC000532	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 10956	CloneNo.: 5F8F1
Source: Mouse	UNIPROT ID: Q13438	
Isotype: IgG1	Full Name: osteosarcoma amplified 9, endoplasmic reticulum associated protein	
Immunogen Catalog Number: AG0106	Calculated MW: 76 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50745-1, OS9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60527-1-PBS. Detection antibody: 60527-2-PBS. Standard:Ag0106. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50745-4, OS9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60527-1-PBS. Detection antibody: 60527-6-PBS. Standard:Ag0106. Range: 0.098-3.125 ng/mL.