

For Research Use Only

OS9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60527-4-PBS



Basic Information

Catalog Number:

60527-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0106

GenBank Accession Number:

BC000532

GeneID (NCBI):

10956

UNIPROT ID:

Q13438

Full Name:

osteosarcoma amplified 9,
endoplasmic reticulum associated
protein

Calculated MW:

76 kDa

Purification Method:

Protein G purification

CloneNo.:

3A11D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50745-2, OS9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60527-3-PBS. Detection antibody: 60527-4-PBS. Standard: Ag0106. Range: 0.195-00 ng/mL



Cytometric bead array standard curve of MP50745-3, OS9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60527-5-PBS. Detection antibody: 60527-4-PBS. Standard: Ag0106. Range: 0.391-100 ng/mL