

For Research Use Only

# OS9 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60527-6-PBS



## Basic Information

Catalog Number:

60527-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0106

GenBank Accession Number:

BC000532

GeneID (NCBI):

10956

UNIPROT ID:

Q13438

Full Name:

osteosarcoma amplified 9,  
endoplasmic reticulum associated  
protein

Calculated MW:

76 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5G12F2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50745-4, OS9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60527-1-PBS. Detection antibody: 60527-6-PBS. Standard: Ag0106. Range: 0.098-3.125 ng/mL