

For Research Use Only

CPT1B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60528-3-PBS



Basic Information

Catalog Number:

60528-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19056

GenBank Accession Number:

BC131570

GeneID (NCBI):

1375

UNIPROT ID:

Q92523

Full Name:

carnitine palmitoyltransferase 1B

(muscle)

Calculated MW:

567 aa, 64 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F12F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50747-2, CPT1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60528-3-PBS. Detection antibody: 60528-4-PBS. Standard: Ag19056. Range: 0.781-100 ng/mL.