

For Research Use Only

# GLRX2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60530-1-PBS



## Basic Information

Catalog Number: 60530-1-PBS	GenBank Accession Number: BC028113	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51022	CloneNo.: 1A10B8
Source: Mouse	UNIPROT ID: Q9NS18	
Isotype: IgG1	Full Name: glutaredoxin 2	
Immunogen Catalog Number: AG34018	Calculated MW: 164 aa, 18 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50749-1, GLRX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60530-1-PBS. Detection antibody: 60530-2-PBS. Standard: Ag34018. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50749-3, GLRX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60530-1-PBS. Detection antibody: 60530-3-PBS. Standard: Ag34018. Range: 0.781-100 ng/mL.