

For Research Use Only

# COX6B2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60531-1-PBS



## Basic Information

<b>Catalog Number:</b> 60531-1-PBS	<b>GenBank Accession Number:</b> BC064548	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 125965	<b>CloneNo.:</b> 1A3B5
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q6YFQ2	
<b>Isotype:</b> IgG1	<b>Full Name:</b> cytochrome c oxidase subunit VIb polypeptide 2 (testis)	
<b>Immunogen Catalog Number:</b> AG33133	<b>Calculated MW:</b> 88 aa, 11 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP50751-1, COX6B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60531-1-PBS. Detection antibody: 60531-2-PBS. Standard:Ag33133. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50751-2, COX6B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60531-1-PBS. Detection antibody: 60531-3-PBS. Standard:Ag33133. Range: 0.781-25 ng/mL



Cytometric bead array standard curve of MP50751-3, COX6B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60531-1-PBS. Detection antibody: 60531-4-PBS. Standard:Ag33133. Range: 0.781-100 ng/mL