

For Research Use Only

GREB1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60532-2-PBS



Basic Information

Catalog Number:

60532-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29941

GenBank Accession Number:

NM_014668

GeneID (NCBI):

9687

UNIPROT ID:

Q4ZG55

Full Name:

GREB1 protein

Calculated MW:

216 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C10G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50752-1, GREB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60532-1-PBS. Detection antibody: 60532-2-PBS. Standard:Ag29941. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50752-3, GREB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60532-4-PBS. Detection antibody: 60532-2-PBS. Standard:Ag29941. Range: 0.098-12.5 ng/mL.