

For Research Use Only

DMWD Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60536-4-PBS



Basic Information

Catalog Number:

60536-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG36845

GenBank Accession Number:

BC019266

GeneID (NCBI):

1762

UNIPROT ID:

Q09019

Full Name:

dystrophia myotonica, WD repeat containing

Calculated MW:

674 aa, 70 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A4A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50758-3, DMWD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60536-4-PBS. Detection antibody: 60536-3-PBS. Standard: Ag36845. Range: 0.098-50 ng/mL