

For Research Use Only

# RTN4IP1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60538-1-PBS



## Basic Information

Catalog Number: 60538-1-PBS	GenBank Accession Number: BC006399	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 84816	CloneNo.: 2F11E4
Source: Mouse	UNIPROT ID: Q8WWV3	
Isotype: IgG1	Full Name: reticulon 4 interacting protein 1	
Immunogen Catalog Number: AG8161	Calculated MW: 396 aa, 44 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50762-1, RTN4IP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60538-1-PBS. Detection antibody: 60538-2-PBS. Standard: Ag8161. Range: 0.098-100 ng/mL.