

For Research Use Only

RTN4IP1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60538-3-PBS



Basic Information

Catalog Number: 60538-3-PBS	GenBank Accession Number: BC006399	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 84816	CloneNo.: 1C4F11
Source: Mouse	UNIPROT ID: Q8WWV3	
Isotype: IgG1	Full Name: reticulon 4 interacting protein 1	
Immunogen Catalog Number: AG8161	Calculated MW: 396 aa, 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50762-2, RTN4P1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60538-3-PBS. Detection antibody: 60538-4-PBS. Standard: Ag8161. Range: 0.195-100 ng/mL.