

For Research Use Only

# RTN4IP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60538-4-PBS



## Basic Information

<b>Catalog Number:</b> 60538-4-PBS	<b>GenBank Accession Number:</b> BC006399	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 84816	<b>CloneNo.:</b> 2B11A3
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q8WWV3	
<b>Isotype:</b> IgG1	<b>Full Name:</b> reticulon 4 interacting protein 1	
<b>Immunogen Catalog Number:</b> AG8161	<b>Calculated MW:</b> 396 aa, 44 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50762-2, RTN4P1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60538-3-PBS. Detection antibody: 60538-4-PBS. Standard: Ag8161. Range: 0.195-100 ng/mL.