

For Research Use Only

FLT3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60540-5-PBS



Basic Information

Catalog Number:

60540-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_004119.

GeneID (NCBI):

2322

UNIPROT ID:

P36888

Full Name:

fms-related tyrosine kinase 3

Calculated MW:

113 kDa

Purification Method:

Protein G purification

CloneNo.:

3C12E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50766-3, FLT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60540-1-PBS. Detection antibody: 60540-5-PBS. Standard:Eg0171. Range: 0.098-100 ng/mL.