

For Research Use Only

PDGFRB Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60542-1-PBS



Basic Information

Catalog Number: 60542-1-PBS	GenBank Accession Number: BC032224	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5159	CloneNo.: 2D4H10
Source: Mouse	UNIPROT ID: P09619	
Isotype: IgG2a	Full Name: platelet-derived growth factor receptor, beta polypeptide	
Immunogen Catalog Number: AG4285	Calculated MW: 1106 aa, 124 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50756-5, PDGFRB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60542-1-PBS. Detection antibody: 60543-1-PBS. Standard: Ag4285. Range: 3.125-200 ng/mL.