

For Research Use Only

PDGFRB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60543-1-PBS



Basic Information

Catalog Number:

60543-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG28376

GenBank Accession Number:

BC032224

GeneID (NCBI):

5159

UNIPROT ID:

P09619

Full Name:

platelet-derived growth factor
receptor, beta polypeptide

Calculated MW:

1106 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

2D5A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50756-5, PDGFRB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60542-1-PBS. Detection antibody: 60543-1-PBS. Standard: Ag4285. Range: 3.125-200 ng/mL.