

For Research Use Only

MRPL27 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60545-1-PBS



Basic Information

Catalog Number: 60545-1-PBS	GenBank Accession Number: BC001066	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51264	CloneNo.: 1A2E5
Source: Mouse	UNIPROT ID: Q9P0M9	
Isotype: IgG1	Full Name: mitochondrial ribosomal protein L27	
Immunogen Catalog Number: AG6870	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50769-1, MRPL27 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60545-1-PBS. Detection antibody: 60545-2-PBS. Standard:Ag6870. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50769-2, MRPL27 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60545-1-PBS. Detection antibody: 60545-3-PBS. Standard:Ag6870. Range: 0.781-100 ng/mL.