

For Research Use Only

# MRPL27 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60545-3-PBS



## Basic Information

Catalog Number: 60545-3-PBS	GenBank Accession Number: BC001066	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51264	CloneNo.: 1C2H10
Source: Mouse	UNIPROT ID: Q9P0M9	
Isotype: IgG1	Full Name: mitochondrial ribosomal protein L27	
Immunogen Catalog Number: AG6870	Calculated MW: 16 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50769-2, MRPL27 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60545-1-PBS. Detection antibody: 60545-3-PBS. Standard: Ag6870. Range: 0.781-100 ng/mL.