

For Research Use Only

SLC39A6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60546-4-PBS



Basic Information

Catalog Number:

60546-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG5676

GenBank Accession Number:

BC039498

GeneID (NCBI):

25800

UNIPROT ID:

Q13433

Full Name:

solute carrier family 39 (zinc transporter), member 6

Calculated MW:

85 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A3D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50771-2, SLC39A6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60546-3-PBS. Detection antibody: 60546-4-PBS. Standard: Ag5676. Range: 0.098-12.5 ng/mL