

For Research Use Only

# CD47 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60551-2-PBS



## Basic Information

Catalog Number:

60551-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC010016

GeneID (NCBI):

961

UNIPROT ID:

Q08722

Full Name:

CD47 molecule

Calculated MW:

323 aa, 35 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H6B12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50813-1, CD47 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60551-2-PBS. Detection antibody: 60551-3-PBS. Standard: Eg31499. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50813-3, CD47 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60551-2-PBS. Detection antibody: 60551-5-PBS. Standard: Eg31499. Range: 0.098-25 ng/mL.