

For Research Use Only

POT1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60556-2-PBS



Basic Information

Catalog Number:

60556-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG36974

GenBank Accession Number:

BC002923

GeneID (NCBI):

25913

UNIPROT ID:

Q9NUX5

Full Name:

POT1 protection of telomeres 1
homolog (S. pombe)

Calculated MW:

71 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H10A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50784-1, POT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60556-1-PBS. Detection antibody: 60556-2-PBS. Standard: Ag36974. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50784-2, POT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60556-2-PBS. Detection antibody: 60556-3-PBS. Standard: Ag36974. Range: 1.563-100 ng/mL.