

For Research Use Only

PTPN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60557-4-PBS



Basic Information

Catalog Number:

60557-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28381

GenBank Accession Number:

BC015660

GeneID (NCBI):

5770

UNIPROT ID:

P18031

Full Name:

protein tyrosine phosphatase, non-receptor type 1

Calculated MW:

435 aa, 50 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B6C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50785-3, PTPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60557-4-PBS. Detection antibody: 60557-5-PBS. Standard: Ag28381. Range: 0.391-100 ng/mL.