

For Research Use Only

PTPN1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60557-5-PBS



Basic Information

Catalog Number: 60557-5-PBS	GenBank Accession Number: BC015660	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5770	CloneNo.: 1G2G2
Source: Mouse	UNIPROT ID: P18031	
Isotype: IgG1	Full Name: protein tyrosine phosphatase, non-receptor type 1	
Immunogen Catalog Number: AG28381	Calculated MW: 435 aa, 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50785-3, PTPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60557-4-PBS. Detection antibody: 60557-5-PBS. Standard: Ag28381. Range: 0.391-100 ng/mL.