

For Research Use Only

SPP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60558-1-PBS



Basic Information

Catalog Number:

60558-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16823

GenBank Accession Number:

BC008938

GeneID (NCBI):

81502

UNIPROT ID:

Q8TCT9

Full Name:

histocompatibility (minor) 13

Calculated MW:

377 aa, 41 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B2B10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50787-1, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-1-PBS. Detection antibody: 60558-2-PBS. Standard:Ag16823. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP50787-2, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-1-PBS. Detection antibody: 60558-3-PBS. Standard:Ag16823. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50787-3, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-1-PBS. Detection antibody: 60558-4-PBS. Standard:Ag16823. Range: 1.563-100 ng/mL