

For Research Use Only

ATP2A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60559-3-PBS



Basic Information

Catalog Number:

60559-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18184

GenBank Accession Number:

BC037354

GeneID (NCBI):

487

UNIPROT ID:

O14983

Full Name:

ATPase, Ca⁺⁺ transporting, cardiac muscle, fast twitch 1

Calculated MW:

1001 aa, 110 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A1E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50788-2, ATP2A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60559-3-PBS. Detection antibody: 60559-4-PBS. Standard: Ag18184. Range: 3.125-100 ng/mL.