

For Research Use Only

RSP01 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60561-2-PBS



Basic Information

Catalog Number: 60561-2-PBS	GenBank Accession Number: BC114966	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 284654	CloneNo.: 2D8D5
Source: Mouse	UNIPROT ID: Q2MKA7	
Isotype: IgG1	Full Name: R-spondin homolog (Xenopus laevis)	
Immunogen Catalog Number: AG34556	Calculated MW: 263 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50792-1, RSPO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60561-1-PBS. Detection antibody: 60561-2-PBS. Standard:Ag34556. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50792-2, RSPO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60561-3-PBS. Detection antibody: 60561-2-PBS. Standard:Ag34556. Range: 0.195-100 ng/mL