

For Research Use Only

LSM14A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60563-4-PBS



Basic Information

| | | |
|---|---|--|
| Catalog Number: 60563-4-PBS | GenBank Accession Number: BC016842 | Purification Method: Protein A Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 26065 | CloneNo.: 3H10D9 |
| Source: Mouse | UNIPROT ID: Q8ND56 | |
| Isotype: IgG2b | Full Name: LSM14A, SCD6 homolog A (S. cerevisiae) | |
| Immunogen Catalog Number: AG12997 | Calculated MW: 463 aa, 51 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50794-2, LSM14A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60563-3-PBS. Detection antibody: 60563-4-PBS. Standard: Ag12997. Range: 6.25-100 ng/mL