

For Research Use Only

RPS8 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60566-1-PBS



Basic Information

Catalog Number: 60566-1-PBS	GenBank Accession Number: BC070875	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6202	CloneNo.: 1B4H2
Source: Mouse	UNIPROT ID: P62241	
Isotype: IgG1	Full Name: ribosomal protein S8	
Immunogen Catalog Number: AG12230	Calculated MW: 208 aa, 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50797-1, RPS8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60566-1-PBS. Detection antibody: 60566-2-PBS. Standard: Ag12230. Range: 6.25-100 ng/mL