

For Research Use Only

RENIN-RECEPTOR,ATP6AP2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60570-1-PBS



Basic Information

Catalog Number:

60570-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17725

GenBank Accession Number:

BC010395

GeneID (NCBI):

10159

UNIPROT ID:

O75787

Full Name:

ATPase, H⁺ transporting, lysosomal accessory protein 2

Calculated MW:

39 kDa

Purification Method:

Protein G purification

CloneNo.:

2B8G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50807-1, RENIN-RECEPTOR,ATP6AP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60570-1-PBS. Detection antibody: 60570-2-PBS. Standard:Ag17725. Range: 1.563-25 ng/mL