

For Research Use Only

COL11A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60571-1-PBS



Basic Information

Catalog Number: 60571-1-PBS	GenBank Accession Number: BC117697	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1301	CloneNo.: 2A4G9
Source: Mouse	UNIPROT ID: P12107	
Isotype: IgG1	Full Name: collagen, type XI, alpha 1	
Immunogen Catalog Number: AG16509	Calculated MW: 1818 aa, 182 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50808-1, COL11A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60571-1-PBS. Detection antibody: 60571-2-PBS. Standard: Ag16509. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50808-2, COL11A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60571-1-PBS. Detection antibody: 60571-3-PBS. Standard: Ag16509. Range: 0.391-100 ng/mL.