

For Research Use Only

# CYP4F2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60573-3-PBS



## Basic Information

Catalog Number:

60573-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33921

GenBank Accession Number:

NM\_001082

GeneID (NCBI):

8529

UNIPROT ID:

P78329

Full Name:

cytochrome P450, family 4, subfamily

F, polypeptide 2

Calculated MW:

60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3F6E2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50811-2, CYP4F2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60573-3-PBS. Detection antibody: 60573-4-PBS. Standard: Ag33921. Range: 0.098-6.25 ng/mL