

For Research Use Only

CYP4F2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60573-4-PBS



Basic Information

Catalog Number: 60573-4-PBS	GenBank Accession Number: NM_001082	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8529	CloneNo.: 3A9C9
Source: Mouse	UNIPROT ID: P78329	
Isotype: IgG3	Full Name: cytochrome P450, family 4, subfamily F, polypeptide 2	
Immunogen Catalog Number: AG33921	Calculated MW: 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50811-2, CYP4F2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60573-3-PBS. Detection antibody: 60573-4-PBS. Standard:Ag33921. Range: 0.098-6.25 ng/mL