

For Research Use Only

CYP2C9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60576-4-PBS



Basic Information

Catalog Number: 60576-4-PBS	GenBank Accession Number: BC020754	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1559	CloneNo.: 1H7G11
Source: Mouse	UNIPROT ID: P11712	
Isotype: IgG1	Full Name: cytochrome P450, family 2, subfamily C, polypeptide 9	
Immunogen Catalog Number: AG9566	Calculated MW: 162aa, 18 kDa; 490aa, 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50815-3, CYP2C9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60576-1-PBS. Detection antibody: 60576-4-PBS. Standard: Ag9566. Range: 0.195-100 ng/mL.