

For Research Use Only

CYP2C9 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60576-5-PBS



Basic Information

Catalog Number: 60576-5-PBS	GenBank Accession Number: BC020754	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1559	CloneNo.: 2B1G7
Source: Mouse	UNIPROT ID: P11712	
Isotype: IgG2a	Full Name: cytochrome P450, family 2, subfamily C, polypeptide 9	
Immunogen Catalog Number: AG9566	Calculated MW: 162aa,18 kDa; 490aa,56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50815-4, CYP2C9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60576-5-PBS. Detection antibody: 60576-2-PBS. Standard: Ag9566. Range: 1.563-100 ng/mL.