

For Research Use Only

UCKL1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60578-4-PBS



Basic Information

Catalog Number: 60578-4-PBS	GenBank Accession Number: BC033078	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 54963	CloneNo.: 1G6C11
Source: Mouse	UNIPROT ID: Q9NWZ5	
Isotype: IgG1	Full Name: uridine-cytidine kinase 1-like 1	
Immunogen Catalog Number: AG10929	Calculated MW: 419 aa, 46 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50820-3, UCKL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60578-4-PBS. Detection antibody: 60578-5-PBS. Standard: Ag10929. Range: 1.563-100 ng/mL.