

For Research Use Only

AMBRA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60582-1-PBS



Basic Information

Catalog Number: 60582-1-PBS	GenBank Accession Number: BC045609	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 55626	CloneNo.: 1F8E11
Source: Mouse	UNIPROT ID: Q9C0C7	
Isotype: IgG1	Full Name: autophagy/beclin-1 regulator 1	
Immunogen Catalog Number: AG4735	Calculated MW: 1208 aa, 133 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50827-1, AMBRA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60582-1-PBS. Detection antibody: 60582-2-PBS. Standard:Ag4735. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50827-2, AMBRA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60582-1-PBS. Detection antibody: 60582-3-PBS. Standard:Ag4735. Range: 0.781-100 ng/mL.