

For Research Use Only

AMBRA1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60582-2-PBS



Basic Information

Catalog Number:

60582-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4735

GenBank Accession Number:

BC045609

GeneID (NCBI):

55626

UNIPROT ID:

Q9C0C7

Full Name:

autophagy/beclin-1 regulator 1

Calculated MW:

1208 aa, 133 kDa

Purification Method:

Protein G purification

CloneNo.:

1B10E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50827-1, AMBRA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60582-1-PBS. Detection antibody: 60582-2-PBS. Standard:Ag4735. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50827-3, AMBRA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60582-2-PBS. Detection antibody: 60582-3-PBS. Standard:Ag4735. Range: 0.781-100 ng/mL.