

For Research Use Only

# PDE3B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60584-3-PBS



## Basic Information

Catalog Number:

60584-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19190

GenBank Accession Number:

BC136565

GeneID (NCBI):

5140

UNIPROT ID:

Q13370

Full Name:

phosphodiesterase 3B, cGMP-  
inhibited

Calculated MW:

1112 aa, 124 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H1C2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50831-2, PDE3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60584-3-PBS. Detection antibody: 60584-4-PBS. Standard: Ag19190. Range: 0.781-50 ng/mL