

For Research Use Only

PTGIS Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60586-2-PBS



Basic Information

Catalog Number: 60586-2-PBS	GenBank Accession Number: BC101809	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5740	CloneNo.: 1C7F2
Source: Mouse	UNIPROT ID: Q16647	
Isotype: IgG1	Full Name: prostaglandin I2 (prostacyclin) synthase	
Immunogen Catalog Number: AG25824	Calculated MW: 500 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50835-1, PTGIS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60586-1-PBS. Detection antibody: 60586-2-PBS. Standard: Ag25824. Range: 0.391-100 ng/mL.