

For Research Use Only

# HSD17B2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60587-3-PBS



## Basic Information

Catalog Number:

60587-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1411

GenBank Accession Number:

BC009581

GeneID (NCBI):

3294

UNIPROT ID:

P37059

Full Name:

hydroxysteroid (17-beta)  
dehydrogenase 2

Calculated MW:

43 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E2B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50836-2, HSD17B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60587-1-PBS. Detection antibody: 60587-3-PBS. Standard:Ag1411. Range: 0.781-100 ng/mL.