

For Research Use Only

# RPN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60589-1-PBS



## Basic Information

|   |  |  |
|---|--|--|
| <b>Catalog Number:</b><br>60589-1-PBS       | <b>GenBank Accession Number:</b><br>BC010839 | <b>Purification Method:</b><br>Protein G Magarose purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>6184                | <b>CloneNo.:</b><br>1B4G5                                      |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>P04843                 |  |
| <b>Isotype:</b><br>IgG1                     | <b>Full Name:</b><br>ribophorin I            |  |
| <b>Immunogen Catalog Number:</b><br>AG28586 | <b>Calculated MW:</b><br>69 kDa              |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP50838-1, RPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60589-1-PBS. Detection antibody: 60589-2-PBS. Standard:Ag28586. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50838-2, RPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60589-1-PBS. Detection antibody: 60589-3-PBS. Standard:Ag28586. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50838-3, RPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60589-1-PBS. Detection antibody: 60589-4-PBS. Standard:Ag28586. Range: 0.098-100 ng/mL