

For Research Use Only

# SERCA2,ATP2A2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60593-1-PBS



## Basic Information

Catalog Number:

60593-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26307

GenBank Accession Number:

BC035588

GeneID (NCBI):

488

UNIPROT ID:

P16615

Full Name:

ATPase, Ca<sup>++</sup> transporting, cardiac  
muscle, slow twitch 2

Calculated MW:

1042 aa, 115 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A5B1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50843-1, SERCA2, ATP2A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60593-1-PBS. Detection antibody: 60593-2-PBS. Standard: Ag5052. Range: 0.098-25 ng/mL.