

For Research Use Only

CACNA1A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60594-3-PBS



Basic Information

Catalog Number:

60594-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG33902

GenBank Accession Number:

NM_000068

GeneID (NCBI):

773

UNIPROT ID:

O00555

Full Name:

calcium channel, voltage-dependent,

P/Q type, alpha 1A subunit

Calculated MW:

282 kDa

Purification Method:

Protein A purification

CloneNo.:

1F7B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50845-2, CACNA1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60594-3-PBS. Detection antibody: 60594-4-PBS. Standard: Ag33902. Range: 3.125-100 ng/mL.